


# medical technologies

surgical supplies



**EkoMed**  
Global

意可曼医疗用品有限公司





# what's inside

4-9	General Surgery
10-15	Cardiovascular
16-21	Cardiology
22-25	Ophthalmology
26-29	Neurosurgery
30-33	Oncology
34-39	Orthopedics
40-47	Gynecology
48-51	E.N.T
52-57	Urology



**EkoMed**  
Global



**EkoMed**  
Global  
意可曼医疗用品有限公司






















 general surgery



UNIVERSAL PACK-A (2 LAYERS) CODE NUMBER:  
GP.110.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	DRAPE WITH ADHESIVE	150 x 240 cm
2	1	DRAPE WITH ADHESIVE	150 x 180 cm
3	2	DRAPE WITH ADHESIVE	75 x 90 cm
4	1	INSTRUMENT TABLE COVER ABSORBENT	150 x 200 cm
5	1	MAYO STAND COVER	78 x 145 cm
6	2	OP TAPE	10 x 50 cm
7	2	SURGDICAL GOWNS REINFORCED	-
8	2	HAND TOWELS	40 x 40 cm

 SMS	 2 layer nonwoven	 SMS	 2 layer nonwoven
 3 layer nonwoven	 Insice film	 3 layer nonwoven	 Insice film
 Blue plastic	 Wrap	 Blue plastic	 Wrap
 Medical adhesive tape	 Transparanty plastic film	 Medical adhesive tape	 Transparanty plastic film



UNIVERSAL PACK-B (2 LAYERS) CODE NUMBER:  
GP.110.002.S

















Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	DRAPE WITH ADHESIVE	150 x 240 cm
2	1	DRAPE WITH ADHESIVE	150 x 180 cm
3	2	DRAPE WITH ADHESIVE	75 x 90 cm
4	1	INSTRUMENT TABLE COVER ABSORBENT	150 x 200 cm
5	1	MAYO STAND COVER	78 x 145 cm
6	2	OP TAPE	10 x 50 cm
8	2	HAND TOWELS	40 x 40 cm






UNIVERSAL PACK-C (3 LAYERS) CODE NUMBER:  
GP.110.003.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	DRAPE WITH ADHESIVE	150 x 240 cm
2	1	DRAPE WITH ADHESIVE	150 x 180 cm
3	2	DRAPE WITH ADHESIVE	75 x 90 cm
4	1	INSTRUMENT TABLE COVER ABSORBENT	150 x 200 cm
	1	MAYO STAND COVER	78 x 145 cm
6	2	OP TAPE	10 x 50 cm
8	2	HAND TOWELS	40 x 40 cm

 SMS	 2 layer nonwoven	 SMS	 2 layer nonwoven
 3 layer nonwoven	 Insicie film	 3 layer nonwoven	 Insicie film
 Blue plastic	 Wrap	 Blue plastic	 Wrap
 Medical adhesive tape	 Transparanty plastic film	 Medical adhesive tape	 Transparanty plastic film



LAPAROTOMY PACK CODE NUMBER:  
GP.112.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	LAPAROTOMY DRAPE	240 x 320 cm
2	1	INSTRUMENT TABLE COVER ABSORBENT	150 x 200 cm
3	1	MAYO STAND COVER	78 x 145 cm
4	2	HAND TOWEL	40 x 40 cm
5	2	OP TAPE	10 x 50 cm























cardiovascular






BY-PASS PACK CODE NUMBER:  
GP.120.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	BY-PASS DRAPE (35 X 42 CM FENESTRATION)	240 x 340 cm
2	1	BY-PASS SPLIT SHEET	200 x 250 cm
3	1	INSTRUMENT TABLE COVER ABSORBENT	150 x 200 cm
4	1	MAYO STAND COVER	78 x 145 cm
5	2	DEFIBRILLATOR BAG	28 x 33 cm
6	2	SURGICAL GOWNS REINFORCED	-
7	2	FOOT COVERS	-
8	4	HAND TOWEL	40 x 40 cm
9	2	OP TAPE	10 x 50 cm

			
SMS	2 layer nonwoven	SMS	2 layer nonwoven
			
3 layer nonwoven	Inside film	3 layer nonwoven	Inside film
			
Blue plastic	Wrap	Blue plastic	Wrap
			
Medical adhesive tape	Transparanty plastic film	Medical adhesive tape	Transparanty plastic film



VALVE REPLACEMENT PACK. CODE NUMBER:  
GP.121.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	VALVE REPLACEMENT DRAPE	240 x 340 cm
2	1	BACK TABLE COVER ABSORBENT	150 x 200 cm
3	1	MAYO STAND COVER	78 x 145 cm
4	2	OP TAPE	10 x 50 cm
5	2	DEFIBRILLATOR BAG	28 x 33 cm
6	2	SURGICAL GOWNS REINFORCED	XL
7	2	TOWEL	40 x 40 cm





THORAX PACK. CODE NUMBER:  
GP.122.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	VALVE REPLACEMENT DRAPE	240 x 340 cm
2	1	BACK TABLE COVER ABSORBENT	150 x 200 cm
3	1	MAYO STAND COVER	78 x 145 cm
4	2	OP TAPE	10 x 50 cm
5	2	DEFIBRILLATOR BAG	28 x 33 cm
6	2	TOWEL	40 x 40 cm



SMS



3 layer nonwoven



Blue plastic



Medical adhesive tape



2 layer nonwoven



Inside film



Wrap



Transparant plastic film









# Cardiology




ANGIOGRAPHY PACK-A CODE NUMBER:  
GP.123.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	FEMORAL ANGIOGRAPHY DRAPE W	200 x 300 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	2	SURGICAL GOWN -STANDARD-	-
4	1	MAYO STAND COVER	78 x 145 cm
5	2	HAND TOWEL	40 x 40 cm
6	1	WRAPPING	100 x 100 cm
7	1	BETADINE SPONGE	-
8	20	GAUZES	7,5 x 7,5 cm



ANGIOGRAPHY PACK-B CODE NUMBER:  
GP.123.002.Sv



Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	RADIAL ANGIOGRAPHY DRAPE W/DUAL TRANSPARENT PANEL	220 x 300 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	2	SURGICAL GOWNS REINFORCED	-
4	1	MAYO STAND COVER	78 x 145 cm
5	2	HAND TOWEL	40 x 40 cm
6	1	WRAPPING	100 x 100 cm
7	1	BETADINE SPONGE	-
8	20	GAUZES	7,5 x 7,5 cm
9	1	SOLUTION BOWL	250 ml
10	1	SOLUTION BOWL	500 ml
11	1	TRAY	2500 ml
12	1	FLUOROSCOPY COVER	R:90 cm
13	1	MANIFOLD KIT	-







ANGIOGRAPHY PACK-C CODE NUMBER:  
GP.123.003.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	ANGIOGRAPHY DRAPE FEMORAL	150 x 300 cm
2	1	INSTRUMENT TABLE COVER PE+ Absorbent	100 x 150 cm
3	2	HAND TOWEL	40 x 40 cm
4	1	FLUOROSCOPY COVER	R:90 cm



ANGIOPLASTY PACK-D CODE NUMBER:  
GP.123.003.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	ANGIOGRAPHY DRAPE FEMORAL	150 x 300 cm
2	1	FLUOROSCOPY COVER	R:90 cm
3	1	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
4	1	SPONGE STICK	-
5	20	GAUZES 8 PLY	
6	2	SOLUTION BOWL	500 ml
7	1	SCALPEL	No: 11
8	2	SURGICAL GOWN REINFORCED	
9	2	HAND TOWEL	40 x 40 cm





# Ophthalmology



**EkoMed**  
Global



















**EkoMed**  
Global  
意可曼医疗用品有限公司






EYE DRAPE PACK-A CODE NUMBER:  
GP.130.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	OPHTHALMOLOGY DRAPE SINGLE LCP	100 x 130 cm
2	1	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
3	2	SURGICAL GOWN	-
4	2	HAND TOWEL	40 x 40 cm


			
SMS	2 layer nonwoven	SMS	2 layer nonwoven
			
3 layer nonwoven	Insice film	3 layer nonwoven	Insice film
			
Blue plastic	Wrap	Blue plastic	Wrap
			
Medical adhesive tape	Transparanty plastic film	Medical adhesive tape	Transparanty plastic film




EYE DRAPE PACK-B CODE NUMBER:  
GP.130.002.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	OPHTHALMOLOGY DRAPE SINGLE LCP	100 x 130 cm
2	1	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
3	2	SURGICAL GOWN	-
4	2	HAND TOWEL	40 x 40 cm





# Neurosurgery





CRANIOTOMY PACK CODE NUMBER:  
GP.140.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	CRANIOTOMY DRAPE WITH LCP	200 x 300 cm
2	4	SIDE ADHESIVE DRAPE	50 x 50 cm
3	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
4	2	HAND TOWELS	40 x 40 cm
5	1	MAYO STAND COVER	78 x1 45 cm
6	2	OP-TAPE	10 x 50 cm



SMS



2 layer nonwoven



3 layer nonwoven



Inside film



Blue plastic




Wrap



Medical adhesive tape




Transparanty plastic film







# Oncology





CHEMOTHERAPY PACK CODE NUMBER:

GP.150.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	TABLE DRAPE Absorbent	80 x 80 cm
2	1	DRAPE	40 x 120 cm
3	1	MASK	-
4	1	PROTECTIVE GLASSES	-
5	1	BONNET	-
6	1-PAIR	CHEMOTHERAPY GLOVES	-
7	1	WASTE POUCH	-
8	2	SURGICAL GOWN REINFORCED 2 PCS	-
9	2	TOWEL	40 x 40 cm



SMS



2 layer nonwoven



3 layer



Inside film



Blue plastic





Wrap




nonwoven  
Medical adhesive tape



Transparently plastic film







# Orthopedics



**EkoMed**  
Global





**EkoMed**  
Global  
意可曼医疗用品有限公司






HIP PACK CODE NUMBER:  
GP.160.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	DRAPE WITH U-SPLIT	200 x 300 cm
2	1	DRAPE WITH ADHESIVE	150 x 240 cm
3	1	DRAPE	75 x 90 cm
4	1	DRAPE	150 x 200 cm
5	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
6	1-PAIR	MAYO STAND COVER	78 x 145 cm
7	1	LEGGINGS	50 x 120 cm
8	1	OP TAPE	10 x 50 cm
9	2	SURGICAL GOWNS REINFORCED	-
10	2	HAND TOWELS	40 x 40 cm




ARTHROSCOPY PACK CODE NUMBER:  
GP.161.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	ARTHROSCOPY DRAPE WITH LIQUID COL LECTION POUCH	200 x 300 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	1	MAYO STAND COVER	78 x 145 cm
4	1	STOCKKINETTE	35 x 85 cm
5	2	SURGICAL GOWN REINFORCED	-
6	4	HAND TOWEL	40 x 40 cm
7	2	OP TAPE	9 x 30 cm












Nr.	Piece	Optimal Extra Contents	Pack Contents Sizes
1	1	CARD BOARD CAMERA COVER	13 x 250 cm
2	1	SKIN MARKER	-
3	1	RULER	-
4	1	ELASTIC BANDAGE	15 x 200 cm
5	1	WRAPPING PAPER	100 x 100 cm



EXTREMITY PACK CODE NUMBER:  
GP.165.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	EXTREMITY DRAPE WITH 7cm ELASTIC APERTURE	200 x 300 cm
2	1	LEGGING	35 x 80 cm
3	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
4	1	MAYO STAND COVER	78 x 145 cm
5	1	HAND TOWEL	40 x 40 cm
6	1	OP TAPE	10 x 50 cm



-  SMS
-  3 layer nonwoven
-  Blue plastic
-  Medical adhesive tape
-  2 layer nonwoven
-  Inside film
-  Wrap
-  Transparant plastic film







Gynaecology



**EkoMed**  
Global




















**EkoMed**  
Global  
意可曼医疗用品有限公司



GYNAECOLOGY PACK CODE NUMBER:  
GP.173.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	LITHOTOMIC GYNAECOLOGY DRAPE WITH LCP	250 x 290 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	2	OP-TAPE	10 x 50 cm
4	2	HAND TOWEL	40 x 40 cm
5	2	SURGICAL GOWN STANDARD	-



 SMS	 2 layer nonwoven	 SMS	 2 layer nonwoven
 3 layer nonwoven	 Inside film	 3 layer nonwoven	 Inside film
 Blue plastic	 Wrap	 Blue plastic	 Wrap
 Medical adhesive tape	 Transparanty plastic film	 Medical adhesive tape	 Transparanty plastic film

DELIVERY PACK CODE NUMBER:  
GP.170.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	SIDE ADHESIVE DRAPE	75x 150 cm
2	1	FLAT BOTTOM DRAPE	75 x 100 cm
3	1	BABY DRAPE	60 x 80 cm
4	2	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
5	2	HAND TOWELS	40 x 40 cm
6	2	REINFORCED SURGICAL GOWNS	L
7	2	LEGGING	





















CAESEREAN PACK-A CODE NUMBER:  
GP.171.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	CAESEREAN DRAPES WITH LCP	240 x 300 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	1	BABY DRAPE	60 x 80 cm
4	2	REINFORCED SURGICAL GOWNS	-
5	2	HAND TOWEL	40 X 40 cm
6	2	OP-TAPE	10 X 50 cm
7	1	UMBILICAL CORD CLAMP	-



 SMS	 2 layer nonwoven	 SMS	 2 layer nonwoven
 3 layer nonwoven	 Inside film	 3 layer nonwoven	 Inside film
 Blue plastic	 Wrap	 Blue plastic	 Wrap
 Medical adhesive tape	 Transparanty plastic film	 Medical adhesive tape	 Transparanty plastic film


CAESEREAN PACK-B CODE NUMBER:  
GP.171.002.S









Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	CAESEREAN DRAPES WITH LCP	150 x 300 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	1	BABY DRAPE	60 x 80 cm
4	2	SURGICAL GOWNS - STANDARD	-
5	2	HAND TOWELS	40 x 40 cm



OPU PACK CODE NUMBER:  
GP.172.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	OPU DRAPE	130 x 150 cm
2	1	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
3	1	PLAIN DRAPE	90 x 90 cm
4	1	TRANSPARENT DRAPE	60 x 60 cm
5	2	LEGGINGS	50 x 120 cm
6	15	GAUZES 12PLY	7,5 x 7,5 cm
7	4	HAND TOWEL	40 x 40 cm



-  SMS
-  3 layer nonwoven
-  Blue plastic
-  Medical adhesive tape
-  2 layer nonwoven
-  Inside film
-  Wrap
-  Transparanty plastic film








E.N.T



**EkoMed**  
Global



**EkoMed**  
Global  
意可曼医疗用品有限公司



E.N.T PACK CODE NUMBER:  
GP.174.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	DRAPE WITH U-SPLIT ADHESIVE 100cm	150 x 240 cm
2	1	FENESTRATED DRAPE WITH SELF-ADHESIVE 10cm. RADIUS	100 x 100 cm
3	2	PLAIN DRAPE	100x 120 cm
4	1	INSTRUMENT TABLE COVER -ABSORBENT	100 x 150 cm
5	2	HAND TOWEL	40 X 40 cm
6	1	WRAPPING PAPER	100 x 100 cm



SMS



2 layer nonwoven



3 layer nonwoven



Insice film



Blue plastic




Wrap




Medical adhesive tape



Transparaty plastic film







# Urology



**EkoMed**  
Global




















**EkoMed**  
Global  
意可曼医疗用品有限公司



TUR PACK CODE NUMBER:  
GP.180.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	T.U.R. DRAPE WITH L.C.P	175 x 240 cm
2	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
3	2	SURGICAL GOWN REINFORCED 2 PCS	-
4	2	HAND TOWEL 2 PCS	40 x 40 cm



			
SMS	2 layer nonwoven	SMS	2 layer nonwoven
			
3 layer nonwoven	Inside film	3 layer nonwoven	Inside film
			
Blue plastic	Wrap	Blue plastic	Wrap
			
Medical adhesive tape	Transparanty plastic film	Medical adhesive tape	Transparanty plastic film

CYSTOSCOPY PACK CODE NUMBER:  
GP.181.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	CYSTOSCOPY DRAPES	75 x 150 cm
2	1	PLAIN DRAPE	75 x 90 cm
3	1	INSTRUMENT TABLE COVER Absorbent	100 x 150 cm
4	1	OP-TAPE	10 x 50 cm
5	2	SURGICAL GOWN REINFORCED	-
6	2	HAND TOWEL	40 x 40 cm
7	2	LEGGINGS	50 x 120 cm





















LAPAROSCOPIC PACK CODE NUMBER:  
GP.182.001.S


Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	LAPAROSCOPY DRAPE WITH LITHOTOMIC LEGGING	175 x 280 cm
2	1	SIDE ADHESIVE DRAPE	75 x 75 cm
3	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
4	2	OP-TAPE	10 x 50 cm
5	2	HAND TOWEL	40 x 40 cm



 SMS	 2 layer nonwoven	 SMS	 2 layer nonwoven
 3 layer nonwoven	 Insice film	 3 layer nonwoven	 Insice film
 Blue plastic	 Wrap	 Blue plastic	 Wrap
 Medical adhesive tape	 Transparanty plastic film	 Medical adhesive tape	 <b>Transparanty</b> plastic film

PERCUTANE PACK CODE NUMBER:  
GP.183.001.S

Nr.	Piece	Pack Contents	Pack Contents Sizes
1	1	PERCUTANE DRAPE WITH L.C.P.	150 x 300 cm
2	1	SCOPY COVER	R: 90 cm
3	1	CAMERA COVER DRAPE	13x 250 cm
4	1	OP-TAPE	10 x 50 cm
5	1	INSTRUMENT TABLE COVER Absorbent	150 x 200 cm
6	2	HAND TOWEL	40 x 40 cm





our Canadian  
headquarters



EkoMed Global, Canada offers its customers many conveniences from storage, sales, finance to door-to-doors shipping. EkoMed Canada Team has more than 20 years of industry experience. They provide "Health Canada" approved products to their Canadian customers.



our sterilization  
facility



In a sterilization facility established on a total area of 2000 square meters in Hangzhou, China, in compliance with advanced technology standards, European Medical Device Directive (MDD) 93/42/ECC certified, in accordance with EO sterilization and ISO class 7-8 clean room standards, sterilization and production of all kinds of medical equipment and medical textile products are carried out at high capacity and speed.



# sterilization

Sterilization comes from the word sterile which is a situation where the biological being loses the capacity to breed or reproduce. Sterilization in its origin is not only limiting the microorganisms but it is also limiting the reproductive capacity of the pathogenic agents which cause infections in the human body.

Any kind of sterilization is in essence realized for not only a point in time but a period of time in our case until the sterilization pouch made-up of polyolefin in some form is a material registered under the trademark of DuPont® as Tyvek™ which has become a generic brand for those medical devices. The content as a medical device is packaged in sterilization pouches in a controlled clean room environment and they are presented for sterilization in various sterilization devices.

There are various types of sterilization listed below:

**A. Steam Sterilization (H<sub>2</sub>O<sub>2</sub> / steam sterilization)**

- Saturated Pure Steam (autoclave)
- Vapor-Phase Hydrogen Peroxide (isolator, clean room)

**B. Dry-Heat Sterilization**

- Sterilization / Depyrogenization Tunnel
- Sterilization / Depyrogenization Oven

**C. Gas Sterilization**

- Ethylene Oxide (autoclave) (flexible chamber gas diffusion ethylene oxide sterilization)
- Formaldehyde (autoclave)
- Gas-Phase Hydrogen Peroxide (isolator)

**D. Ionizing Radiation**

- GAMA (Cobalt 60, radio isotope decay)
- X-ray Radiation
- Electron beam radiation

**E. Filtration**

- Sterilization Grade Filter (0.2 um filter- liquid, gas / air)



Our type of ethylene oxide sterilization has been evaluated as a new generation of Sterilizers known as Flexible Chamber Ethylene Oxide Gas Diffusion System sterilizers and is a high-tech sterilizer that completes the sterilization process by using low temperature and low dose of Ethylene Oxide gas. The technology has eliminated the infrastructure requirements such as compressed air, water, vacuum lines and drainage channels required by conventional gas sterilizers. Cartridges/ ampoules containing Ethylene Oxide gas are used for sterilization.





**EkoMed**  
Global  
意可曼医疗用品有限公司



Canada Headquarters

5510, Mainway, Burlington, Ontario  
Canada, L7L6C4  
Telephone: +1(289)-846-3444  
[sales@ekomedglobal.com](mailto:sales@ekomedglobal.com)  
[www.ekomedglobal.com](http://www.ekomedglobal.com)



China Manufacturing

TONGLU / HANGZHOU



[www.ekomedglobal.com](http://www.ekomedglobal.com)

